



Carterra® LSA™ Instrument

Disruptive Technology for mAb
Screening and Characterization



Fully Integrated High Throughput
SPR Platform


Disruptive Technology for mAb Screening and Characterization

Understanding a monoclonal antibody's (mAb) mechanism of action (MOA) is fundamental to the discovery of superior therapeutics because a mAb's epitope largely dictates its biological function.

The LSA is the only fully integrated, high-throughput mAb screening and characterization platform. The system combines patented flow printing microfluidics with surface plasmon resonance (SPR) detection to deliver high-throughput kinetic and epitope analysis workflows that support state-of-the-art mAb discovery programs.

Automated flow cell switching between Single Flow Cell and 96-Channel Printhead

- Unique fluidics integrates both single flow cell and 96-Channel Printhead modes with high throughput SPR
- Up to 384 reaction spots + 48 reference spots per array
- Supports capture formats and standard amine coupling



Deliver analyte to as many as 384 ligands immobilized on a single array in a minimum volume using the Single Flow Cell.

Immobilize up to 384 ligands on a single array using patented Flow Printing technology.

Unrivaed Throughput for Key mAb Discovery Applications

- **Kinetics & Affinity:**
 - » **Capture kinetics:** Screen up to 1,152 mAbs in a single assay
 - » **Coupled kinetics:** Up to 384 immobilized ligands analyzed simultaneously
- **Epitope Binning:**
 - » **Interrogate up to 384×384 mAbs in a single assay**


- **100x the data**
- **In 10% of the time**
- **With 1% of the sample requirements of other systems**

Highest Throughput SPR Results Without Compromising on Data Quality

Other leading systems sacrifice throughput for data quality, but you no longer have to restrict your R&D based on technology limitations.

The LSA not only provides the highest throughput of any SPR system on the market, it does so without any loss in data quality.

- Excellent agreement in kinetic rate constants
- Data correlates across wide affinity range <100pM to >100nM
- LSA consumes 1% sample of Biacore™ 8K
- LSA analyzes 384 binding interactions in a single day/run
- LSA has powerful batch-mode fitting software to facilitate analysis



Powerful, Intuitive Software Integrates mAb Discovery

The streamlined **Navigator User Interface** enables quick and efficient experimental setup for a diverse range of experiments using intuitive applications for each of the core mAb discovery areas of interest, thereby minimizing hands-on time. Dedicated Kinetics and Epitope data analysis packages provide rapid evaluation, analysis, and visualization of large data sets; we use patented software tools for multiple data views that aid the discovery of unique high-value mAbs.




The **Kinetics** data analysis platform is built to rapidly handle 1000's of interactions in a matter of seconds, using batch-processing routines to simplify the process and speed access to final fitted data. The software automatically applies QC to flag poorly-performing clones and facilitates multiplexed studies of mAbs targeting multiple distinct antigens. Kinetic data can be viewed as 384 tile plots, or iso-affinity plots, with each individual curve and raw data only a click away.



Simultaneous kinetic analysis of 384 antigen/antibody binding interactions using high throughput SPR (Left), with (Right) Enlarged view of the data from three spots showing antibodies that bound their target with diverse kinetics (low, medium, and high affinities, from top to bottom).

The **Epitope** software enables rapid and efficient binning of up to 384x384 mAb competition matrix experiments that reveal exquisite epitope differentiation and identify unique/nuanced binders. Data is presented across three visualization panels to provide a comprehensive view of raw data, sorted heat maps and easy-to-interpret network plots, as shown below.



Training

The LSA is supported by a variety of training options that provide customers with the flexibility to choose the level most appropriate to their needs:

Basic New User Training

- Introduction to the Carterra instrument, with hands-on training covering the basic principles of label-free SPR analysis, applications, and general system management and maintenance. (2 days, On-Site)

Advanced User Training

- Advanced overview of the Carterra instrument, including pre-training consultation, design and development of custom assays, hands-on training, data analysis, and post-training consultation. (2 days, On-Site)

Refresher Basic Training

- Introduction to the Carterra instrument that covers applications, assay development, and provides an overview of data analysis. Via online platform or in-person at Carterra's Customer Experience Center. (1 day, Remote)

Advanced Applications Training

- Review of the Carterra platform including an introduction to the system, reviews that applications and workflows, and provides in-depth training on advanced assay development, assay support, and data analysis. This course includes a pre-training and post-training consultation. (4 days, On-Site)

Kinetics Training

- Introduction to the Kinetics application, including assay development, optimization, and data analysis. (1 day, On-site or Remote)

Epitope Binning Training

- Introduction to the binning function of the Epitope application, including assay development, optimization, and data analysis. (1 day, On-site or Remote)

Custom Applications and Training

- Custom applications development and training defined by customer. (1 day On-site or ½ day Remote)

Sensor Chips

Visit our [online store](#) to get more information about our sensor chips, consumable offerings and to place orders.

Name	Description
HC30M	Polycarboxylate hydrogel, medium charge density, 30nm coating thickness
HC200M	Polycarboxylate hydrogel, medium charge density, 200nm coating thickness
CMDP	2D planar carboxymethyl-dextran surface <5nm coating thickness
CMD50M	Carboxymethyl-dextran hydrogel, medium charge density, 50nm coating thickness
CMD200M	Carboxymethyl-dextran hydrogel, medium charge density, 200nm coating thickness
SAD200M	Streptavidin, immobilized in a carboxymethyl-dextran hydrogel, medium charge density, 200nm coating thickness
NiHC200M	Poly - NTA derivatized linear polycarboxylate hydrogel, medium charge density, 200nm coating thickness
PAHC200M	Protein A derivatized linear polycarboxylate hydrogel, medium charge density, 200nm coating thickness
PAGHC30M	Protein A/G derivatized linear polycarboxylate hydrogel, medium charge density, 30nm coating thickness
PAGHC200M	Protein A/G derivatized linear polycarboxylate hydrogel, medium charge density, 200nm coating thickness

