

Rapid Kinetic and Antagonistic Profile Screening of Hybridoma Supernatants, a COVID-19 Case Study

This Application Note was developed in collaboration with Takis Biotech, Rome, Italy.

Key Takeaways

- Perform high-throughput capture kinetics and competition characterization on hybridoma samples.
- Screen 1,800 clones using only 70 µg of sample.
- Analyze broad kinetic ranges from a single concentration series of samples.
- Quickly view and select the most suitable clones with user-friendly analysis software.

Introduction

The rapid global spread of infectious agents necessitates the urgent need to identify therapeutic candidates quickly. However, the modern drug discovery process is challenged by time constraints from the analysis, characterization, selection and development of lead therapeutic candidates.

Here we present the rapid kinetic and functional screening of 1,800 hybridoma cell culture supernatant samples containing SARS-Cov-2 receptor-binding domain (RBD) specific antibodies generated from immunized mice. Following the full kinetic characterization of the antibodies, a selection of high affinity antibodies antagonizing the interaction between the RBD and Angiotensin-Converting Enzyme 2 (ACE2), which is required for viral entry¹, has been performed using the Carterra LSA platform. The whole process including data evaluation, analysis and clone selection has been performed within a work week highlighting the benefits implementing high-throughput SPR (HT-SPR) in early drug discovery. Carterra's LSA platform with proprietary Kinetics and Epitope software is enabling drug discovery and contract research organizations (CROs) to streamline the discovery process by generating deep

information content data. This is resulting in making decisions earlier in the discovery process thereby facilitating stronger selection of therapeutic antibody candidates to further investigate.

This COVID-19 study has been performed in collaboration with Takis Biotech whom we thank for supplying the 1,800 hybridoma culture supernatant containing potentially COVID-19 neutralizing antibodies.

Method

Surface Preparation of the Capture Lawn (Anti Mouse-Fc mAb, Capture Reagent on the LSA)

Antibodies were captured via their FC domain using immobilized anti mouse-Fc mAb, which was coated as a lawn onto an HC200M sensor chip (Carterra, p/n 4287) using standard EDC/sulfo-NHS-mediated amine coupling method in a running buffer of HBSTE (10 mM Hepes pH7.4, 150 mM NaCl, 3 mM EDTA, 0.01% Tween-20 (Carterra p/n 3631). The LSA single flow cell was used to prepare the lawn, which involved (1) activating the chip with a 10 minutes injection of a mixture of 0.133 M EDC + 0.033 M sulfo-NHS in 0.1 M MES pH 5.5 (Carterra p/n 3626), (2) coupling 75 µg/ml goat anti mouse-Fc mAb (Jackson ImmunoResearch

#115-005-071) in 10 mM sodium acetate pH 4.5 (Carterra p/n 3802) for 12 minutes and (3) blocking excess reactive esters with a three minutes injection of 1 M ethanolamine HCl pH 8.5 (Carterra p/n 3632). This yielded final coupled levels of about 5,000 RU.

Surface Preparation of the Array (Hybridoma Supernatants)

1,800 hybridoma cell culture supernatants were supplied in 96-well plates. Samples for capture were prepared in 384-deep-well plate by performing a five-fold dilution using HBSTE running buffer (final volume 220 µL). A total of five 384-well plates were prepared. Samples were flow printed for 15 minutes in batches of 96 supernatants in parallel using the 96-channel printhead, which passed a fixed volume of sample (200 µL) back-and-forth across the chip surface, thereby enriching the capture of antibodies onto discrete spots. An entire 384-array of spotted antibodies was prepared by serially docking the printhead onto each of the four nested print group locations. At the end of the antibody capture step, samples are flown back in their well of origin authorizing to re-use the same sample solutions to repeat the measurement or to perform further kinetic characterization with additional antigens.

Capture Kinetics

The target antigen RBD was supplied as a purified recombinant monomer and prepared in running buffer of HBSTE + 0.5 mg/ml BSA as an eight points three-fold dilution series with concentration starting at 1 μ M and final volume of 300 μ L.

Analytes were delivered in the single flow cell in ascending concentration, allowing a 5 minutes association time and a 15 minutes dissociation time. Six binding cycles of buffer alone were injected before the antigens samples to stabilize the surface and provide blanks for double-referencing the data.

The last RBD antigen injection (1 μ M) was followed by an ACE2 protein injection at 500 nM in order to determine the antagonistic profile of each antibody. Binding data were analyzed in Carterra's Kinetics data analysis software to determine kinetic rate and affinity constants by fitting a simple Langmuir model to the double-referenced data.

Results

Chip Type

For screening purposes a high capacity chip (HC200M) was used. High surface density combined with bi-directional sample delivery across spots on a chip surface allows for the enrichment of capture levels, even from low titer antibody samples from low expressing clones.

Non-Regenerative Kinetics

The simultaneous kinetic characterization of 384 antibodies required 14 μ g of RBD antigen. The measurement of the whole panel of 1,800 samples was performed using 70 μ g of antigen as well as 44 μ L of each culture supernatant.

A 384-mAb panel using 460 pM to 1 μ M monomeric antigen as analyte, is shown as a tile view in **Figure 1**. Each panel represents the binding responses (colored by analyte concentration, with a blue/green palette) and global fit (in red) obtained for the antigen interacting with a single mAb. Ligands with low or no binding were automatically flagged by the software and are shown as grey and the ones displaying deviation from the fit model are shown in yellow.




Figure 1. Capture kinetics was performed on 384 ligands simultaneously using the LSA. The specific monomeric antigen was injected as analyte from 0.46 nM to 1 μ M over the entire ligand set.

Figure 4. (Left) Assessment of the ability of antibodies to inhibit the binding of ACE2 to RBD.




Figure 5. Signals recorded for a RBD titration (left) followed by a 500 nM ACE2 injection (right). Antibody A displays a high affinity for RBD but does not block the binding of ACE2 to RBD. Antibody B displays a high affinity for RBD and blocks the binding of ACE2 to RBD.

Summary

In the present study, we were able to characterize the kinetics of 1,800 samples in 48 hours with 70 μ g of antigen. Additionally, evaluation of large data sets and performing clone selection is quick and easy using Carterra's powerful and intuitive Kinetics software. The whole screening campaign, including assay measurement and data evaluation, was completed within a week. The Carterra LSA enables the kinetic characterization of biomolecular interactions in a 384-array format and is well suited to be integrated into an early discovery workflow. Non-regenerative capture kinetic experiments performed by HT-SPR take significantly less time and consume dramatically less antigen than conventional biosensor approaches.

References

1. Tai W, He L, Zhang X, Pu J, Voronin D, Jiang S, et al. Characterization of the receptor-binding domain (RBD) of 2019 novel coronavirus: implication for development of RBD protein as a viral attachment inhibitor and vaccine. *Cell Mol Immunol*. 2020, 17, 613
2. Shang J, Ye G, Shi K, Wan Y, Luo C, Aihara H, et al. Structural basis of receptor recognition by SARS-CoV-2. *Nature* 2020, 581, 221

Carterra technology is protected by the following patents and other patents pending:
8,210,119, 8,211,382, 8,383,059, 8,999,726, 9,682,372, 9,682,396, 10,825,548

© 2021 Carterra Inc. Carterra, the Carterra Logo and all other trademarks are property of Carterra unless otherwise specified.

AN109-REV012