



Software Version:
1.8.0.3470

Release Date:
3.28.2022

Contents

- Changes Since 1.7.1.3055 Release
 - Features
 - Updates
 - Known Issues
 - Installation instructions

1 Changes Since Last Release

1.1 Features

- Success dialog for Save and Export
- Improved functionality of error reporting dialog

1.2 Updates

- The headers in the Ligands and Analytes tables aren't displayed in full Peptide Mapping mode
- The new experiment retains the navigation from the previous experiment in Peptide and Mutant Mapping
- An unstable error appears when opening Premix Binning .ebp file
- Items to export aren't displayed in Peptide Mapping mode
- Competition thresholds are not properly sorting competition
- Competition threshold reverts to previous value





1.3 Known Issues

- An unstable error appears when opening a .ebp files (originally .Isadata file) in Peptide Mapping mode
- Community clustering can force convergence of dissimilar epitope bins

2 Installation/Upgrade Instructions

- i. Run the Epitope installer executable file
- ii. Accept User Account Control permissions
- iii. Complete standard installation steps as instructed
- iv. Open application from Windows menu or desktop shortcut

