



# Software Version: 1.8.1.3704

Release Date: 8.22.2022

---

## Contents

- Changes Since 1.8.0.3470 Release
  - Features
  - Updates
  - Known Issues
  - Installation instructions

## 1. Changes Since Last Release

### 1.1 Features

- Release Notes in File menu now link to Carterra Customer Portal

### 1.2 Updates

- Resolved error with application closing while navigating to the Heat Map

### 1.3 Known Issues

- An unstable error appears when opening a .ebp files (originally .lsadata file) in Peptide Mapping mode
- Community clustering can force convergence of dissimilar epitope bins

## 2. Installation/Upgrade Instructions

- i. Run the Epitope installer executable file
- ii. Accept User Account Control permissions
- iii. Complete standard installation steps as instructed
- iv. Open application from Windows menu or desktop shortcut

THIS DOCUMENT MAY CONTAIN CONFIDENTIAL AND PROPRIETARY INFORMATION.  
ANY UNAUTHORIZED REVIEW, USE, DISCLOSURE OR DISTRIBUTION IS PROHIBITED

Carterra technology is protected by the following patents and other patents pending:  
8,210,119, 8,211,382, 8,383,059, 8,999,726, 9,682,372, 9,682,396, 10,825,548

© 2022 Carterra Inc. Carterra, the Carterra Logo and all other trademarks are property of Carterra unless otherwise specified.

RN1.8.1.3704-REV002



info@carterra-bio.com

Headquarters:  
825 N 300 W Suite C309  
Salt Lake City, UT 84103

[www.carterra-bio.com](http://www.carterra-bio.com)