

Recommended Consumables for Carterra[®] LSA

The following consumables are recommended to achieve optimal performance on your LSA instrument. These consumables are not available directly from Carterra.

	Vendor	Type	Part Number	Max Volume
Microtiter Plates	Greiner	96 standard	650201	300 µL
	Greiner	96 deep	780201	1100 µL
	Greiner	384 deep	781270	225 µL
	Cole-Parmer	384 extra deep	EW-07904-33	300 µL
Vials	Simport or VWR	Microtubes	T341-4TTP	1500 µL
	Simport or VWR	Caps with septa	T347AQX	N/A
Conical Tubes	Syringa	15 ml with septa	1J711 [‡]	15 ml
	Syringa	50 ml with septa	1J911 [‡]	50 ml
Reagents	Thermo Fisher Scientific	Sulfo NHS	24510	N/A
	Thermo Fisher Scientific	Premium grade EDC	PG82079	N/A
Glass Slides	Electron Microscopy Sciences	1.2 mm thick	71880-01	N/A

[‡] These items may not be available in all countries or regions. Sapidyne is the official distributor for Syringa and may carry a comparable part.

Disclaimer: Carterra strongly recommends using these consumables with the LSA instrument to achieve optimal performance. Carterra makes no claims to the pricing or availability of these vendor consumables.

Please direct any inquiries to questions@carterra-bio.com